(12) PATENT APPLICATION PUBLICATION

(19) INDIA

(22) Date of filing of Application :15/03/2024

(43) Publication Date : 19/04/2024

(54) Title of the invention : A DENGUE REVERSE TRANSCRIPTION POLYMERASE CHAIN REAC TION (RT-PCR) ASSAY KIT FOR DETECTION AND SEROTYPING OF DENGUE VIRUS

(51) International classification	:C12Q0001700000, C12Q0001685100, C12Q0001688600, C12Q0001689000, C12Q0001686000	(71)Name of Applicant :1)JSS Academy of Higher Education & Research (Deemed to be university
(86) International Application No Filing Date	:NA :NA	Address of Applicant :Sri Shivarathreeshwara Nagara, Mysuru, Karnataka, India-570015 Mysore Name of Applicant : NA
(87) International Publication No	: NA	Address of Applicant : NA (72)Name of Inventor :
(61) Patent of Additio to Application Number Filing Date	n:NA er :NA	1)Dr. Sumana M N Address of Applicant :JSS Academy of Higher Education and Research, Mysore, Karnataka, India Mysore
(62) Divisional to Application Number Filing Date	:NA :NA	2)Ravindran Kumar Address of Applicant :JSS Academy of Higher Education and Research, Mysore, Karnataka, India Mysore

(57) Abstract :

A DENGUE REVERSE TRANSCRIPTION POLYMERASE CHAIN REAC TION (RT-PCR) ASSAY KIT FOR DETECTION AND SEROTYPING OF DENGUE VIRUS ABSTRACT 5 A dengue reverse transcription polymerase chain reaction (RT-PCR) assay kit for de tection and serotyping of dengue virus is provided. The dengue RT-PCR assay kit in cludes a plurality of primer sets for dengue virus detection and serotyping including at least one forward primer and reverse primer, a synthetic positive control sequence for 10 validating the effectiveness of the RT-PCR assay kit, an internal control sequence for monitoring an efficiency of an entire process and identifying potential issues during sample processing and amplification, one or more standards for serving as benchmarks for accurate quantification and spanning copy numbers from 10^5 to 10^1 and an Eva green dye and high-resolution melt analysis (HRMA) for cost-effective and sensitive 15 detection. The Eva green dye binds to double-stranded DNA, facilitating realtime monitoring of amplification. [FIG.1]

No. of Pages : 33 No. of Claims : 14